

Product Datasheet

PEP3 Antibody, pAb, Rabbit**RBP1303, RBP1304****LUCERNA-CHEM**

www.lucerna-chem.ch

Catalog No.	RBP1303, RBP1304
Product name	Anti-PEP3 antibody
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Amount/Conc.	100ug/100ul
Purification	Immunoaffinity chromatography
Conjugation	unconjugated
Specificity	UniProt / Swiss-Prot No.: Q96EH5, 60S ribosomal protein L39-like PEP3 AA sequence: WIQMKPGSKIRYN MSSHKTFTIKRFLAKKQKQNRPIPQWIQMKPGSKIRYNSKRRHWRRTKLG
Immunogen	The antiserum was produced against synthesized peptide derived from customers Sequence: -WIQMKPGSKIRYN
Antigen species	?
Reactivity	?
Crossreactivity	?
Tested applications	ELISA
Other applications	Working concentrations for specific applications should be determined by the investigator. The appropriate concentrations may depend on secondary antibody affinity, antigen concentration, detection method, temperature, incubation length, and other factors. The suitability of this antibody for applications other than those listed above has not been determined.
Storage Buffer	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
Stability	Stability at -20°C / 1 year
Regulation Status	For research use only (RUO)