


HRP-Conjugate/Protein Stabilizer

AS15 HRPstab | Bovine and mammalian protein free HRP-conjugate and protein stabilizer (10 ml, 100 ml, 1 L)

One-component, liquid, ready-to-use formula

Agrisera HRP-conjugate stabilizer preserves and stabilizes the activity of HRP conjugates, proteins and other biomolecules. It is designed as a ready-to-use stabilizer, buffer, diluent and possible coat block for immunoassays. This stabilizer does not contain mammalian proteins, like bovine proteins of BSA. This minimizes background issues, unwanted protein interactions and concerns with possible transmission of bovine spongiform encephalopathy (BSE).



Advantages

- ✓ No mammalian proteins included
- ✓ Minimizes background interactions
- ✓ Blocking and stabilization at once
- ✓ Can be used directly as an assay buffer
- ✓ Stable for at least 2 years
- ✓ Stable through freeze cycles and at RT
- ✓ Not affected by glycerol
- ✓ Produced based on SureGuard chemical technology

Applications

- ✓ ELISA
- ✓ Western blot
- ✓ Southern blot
- ✓ Northern blot
- ✓ Lateral flow strips

Related products

Product name.....	Product number
AgriseraELISA SuperBright (10 ml, 100 ml, 1 L).....	AS18 ELISA-ECL
AgriCoat-Block Non mammalian Protein Coat Block (5X) (10 ml, 100 ml, 1 L).....	AS18 AgriCoat-Block
AgriseraTMB Diluent (100 ml).....	AS15 Diluent
AgriseraTMB HRP Substrate (10 ml, 100 ml, 1 L).....	AS15 TMB-HRP