






Bioss
ANTIBODIES

AbBy Fluor[®] Dyes




Bright and stable fluorochromes




Fluorochrome guide

	Common Filters (nm)	Bioss Dyes	Other Dyes			
 Ultraviolet 350 nm	450/50	AbBy Fluor® 350	Alexa Fluor 350			
	515/30		BD Horizon BUV496	Qdot 525		
	530/43		Qdot 545			
	585/42		Qdot 565	Qdot 585		
	610/20		Qdot 605			
	630/69		Qdot 625			
	660/20		Qdot 655			
	695/40		Qdot 705			
	780/60		Qdot 800			
 Violet 405 nm	450/50	AbBy Fluor® 405	Alexa Fluor 405	BD Horizon BV421	BD Horizon BV480	BD Horizon V450
			Calcein Violet 450 AM	CoraLite 405S	DyLight 405	eBFP
			eCFP	eFluor 450	iFluor 405	PromoFluor 405
			StarBright Violet 440	StarBright Violet 475	SuperBright 436	SuperBright 780
			SuperNova v428	ViaFluor 405	Vio Bright V423	VioBlue
			violetFluor 450	Zombie Violet		
	510/50		BD Horizon BV510			
	525/50		BD Horizon V500	BD Horizon BV510	VioGreen	violetFluor 500
			eFluor 506	Spark Violet 538	StarBright Violet 515	
	530/43		Alexa Fluor 430			
	550/40		Cascade Yellow			
	585/40		StarBright Violet 570	Zombie Yellow		
	585/42		BD Horizon BV570	Pacific Orange		
	610/20		BD Horizon BV605	StarBright Violet 610	SuperBright 600	SuperNova v605
	660/20		BD Horizon BV650	StarBright Violet 670	SuperBright 645	
695/40		StarBright Violet 710				
710/50		BD Horizon BV711	SuperBright 702			
780/60		BD Horizon BV750	BD Horizon BV785	BD Horizon BV786	StarBright Violet 760	
			StarBright Violet 790	SuperNova v786		
 Blue 488 nm	510/50		GFP			
	515/30		AmCyan	CoraLite 488	SureLight 488	
	525/30		ViaFluor 488			
	525/50		Calcium Green	DyLight 488	eGFP	iFluor 488
			NovaFluor Blue 510	NovaFluor Blue 530	Vio515	VioBright 515
			VioBright FITC	YFP		
	530/43	AbBy Fluor® 488 FITC	Alexa Fluor 488	Alexa Fluor 514		
	575/36		NovaFluor Blue 555			
	585/42		NovaFluor Blue 585	Spark Blue 574		
	610/20		NovaFluor Blue 610 / 30S		NovaFluor Blue 610 / 70S	
	660/20		NovaFluor Blue 660 / 120S		NovaFluor Blue 660 / 40S	
	670/30	PerCP	PE-Vio770	PerCP-Vio700		
695/40	PerCP-Cy5.5	PerCP-eFluor 710	StarBright Blue 700	SureLight PerCP Cyanine 5.5		
 Yellow/ green 532 nm	550/40		iFluor 514	iFluor 532		
	575/36	AbBy Fluor® 555	Alexa Fluor 532	Alexa Fluor 546	Alexa Fluor 555	CoraLite 532
	585/42		eFluor 565NC	Spark YG 570		
	586/14		Cy3			

Fluorochrome guide

	Common Filters (nm)	Bioss Dyes	Other Dyes			
	575/36		CoraLite 555			
	582/15		iFluor 555	PromoFluor 555P		
	585/42		CoraLite 568	Cy3.5	Cyanine 555	DyLight 549
			DyLight 550	Dyomics 547	eFluor 570	eFluor 585NC
			MaxLight 550	NovaFluor Yellow 590	PE-Dyomics 590	
	586/14	PE	SureLight RPE			
	586/25		Janelia Fluor 549	NovaFluor Yellow 570	Oyster 550	
	610/20		eFluor 605NC	Janelia Fluor 585	NovaFluor Yellow 610	PE-CF594
			PE-Dazzle 594	PE-DyLight 594	PE-eFluor 610	PE-Vio615
	630/69		Alexa Fluor 568	PE-Alexa 610	PE-ATTO 594	PE-Fire 640
			PE-Texas Red			
660/20		NovaFluor Yellow 660	PE-Alexa 647			
670/30	PE-Cy5	PE-Cy5				
695/40	PE-Cy5.5	NovaFluor Yellow 690	NovaFluor Yellow 700	PE-Alexa 700	PE-Cy5.5	
		PE-Dyomics 647	PE-Fire 700	SureLight PE Cyanine 5.5		
780/60	PE-Cy7	PE-Alexa 750	PE-Cy7	PE-Fire 780	SureLight PE Cyanine 7.0	
	610/20		CoraLite +595	CoraLite 594	eFluor 615	Flash Phalloidin Red 594
			iFluor 594	PhenoVue Fluor 594	PromoFluor 594	Texas Red
	630/69	AbBy Fluor® 594	Alexa Fluor 594	Alexa Fluor 610	Alexa Fluor 633	
	660/20	AbBy Fluor® 647	Alexa Fluor 647	CoraLite 640	CoraLite 647	CoraLite 647A
		APC	DyLight 649	Dyomics 647	iFluor 633	iFluor 647
		Cy5	IRDye 650	NovaFluor Red 660	SureLight APC	
	695/40	AbBy Fluor® 680	Alexa Fluor 660	Alexa Fluor 680	APC-Alexa 680	APC-Alexa 700
		APC-Cy5.5	APC-R700	CoraLite +650	CoraLite 680R	DyLight 650
		Cy5.5	DyLight 680	iFluor 680	iFluor 700	NovaFluor Red 685
			NovaFluor Red 700	NovaFluor Red 710	ViaFluor 647	
725/40		Vio Bright R720	Vio R720			
730/45		Alexa Fluor 700				
780/60	AbBy Fluor® 750	Alexa Fluor 750	APC-Alexa 750	APC-CF750	APC-CF750T	
	APC-Cy7	APC-CF770	APC-eFluor 780	APC-Fire 750	APC-Fire 810	
		APC-H7	APC-Vio770	Cyanine 7 SE		

We are excited about the launch of AbBy Fluor® dyes. These bright and stable fluorochromes are supplied under Bioss brand. Their participation in support of multiplex immunoassays will dramatically facilitate your panel design.



Bioss Inc. | 300 Tradecenter Drive, Suite 4610, Woburn, MA 01801 | +1.781.569.5821 | support@biossusa.com