

ELISA Kits & Assay Development

ELISA Kits & Assay Development

ICL R&D to IVD

What is ELISA?

Enzyme-linked immunosorbent assay (ELISA) is one of the most widely used immunoassay formats due to its sensitivity, flexibility, and ability to provide analyte quantitation. Additionally, the relatively low cost of ELISA and its adaptability to high-throughput testing make it the ideal platform for many different research and in vitro diagnostic (IVD) needs.

ICL's ELISA kits are designed to help detect a broad range of human and veterinary targets in sample types ranging from tissue culture supernatants to milk. Plus, they are backed by the support of our highly experienced team and 40+ years of industry knowledge.

Unparalleled Reliability

At ICL, we offer a broad range of ready-to-use ELISA kits and custom ELISA development for human and animal research, ensuring reliable research results every time. We also offer an extensive selection of antibodies and protein reference standards for researchers developing ELISA in-house. And because our products are manufactured in the USA, fast and efficient delivery is always available.



Assay development reagents for use in combination with our Antibodies and Antigens

Antibodies for LFD - Positive Control
Anti-BSA Antibodies
Anti-Chicken IgY h+l Antibodies
Anti-Goat IgG h+l Antibodies
Anti-Human IgG h+l Antibodies
Anti-Mouse IgG h+l Antibodies
Anti-Rabbit IgG h+l Antibodies
Anti-Rat IgG h+l Antibodies

Purified IgG/IgY
Canine IgG h+l
Chicken IgY h+l
Feline IgG h+l
Goat IgG h+l
Mouse IgG h+l
Rabbit IgG h+l
Rat IgG h+l

Normal Sera
Normal Canine Sera
Normal Chicken Sera
Normal Feline Sera
Normal Mouse Sera
Normal Rabbit Sera
Normal Rat Sera

Miscellaneous Antibodies
Anti-Protein A Antibodies
Anti-Protein G Antibodies
Anti-Protein L Antibodies
Anti-FITC Antibodies
Anti-Alex Fluor 647 Antibodies
Anti-6xHIS Tag Antibodies
Anti-C-MYC Tag Antibodies
Anti-Green Fluorescent Protein Antibodies
Anti-GST Tag Antibodies
Anti-HA Tag Antibodies

Featured Products

ICL R&D to IVD

Preformatted ELISA kits with specificity to Human and Veterinary targets	
Albumin	IgG1, IgG2A, IgG2B, IgG2C, IgG3
Alpha 1-Acid Glycoprotein (AGP)	Kappa
Alpha 1-Antichymotrypsin	KIM-1
Alpha 1-Antitrypsin	Lactoferrin
Alpha 2-Macroglobulin	Lambda
Antithrombin III	L-FABP
ApoH (beta 2-Glycoprotein)	LRG1
Apolipoprotein A1	Myoglobin
Beta 2-Microglobulin	NGAL
C1q	Osteopontin
C3	Ovotransferrin (Conalbumin)
Canine Soluble IL-2 Receptor (sIL2-R)	Plasminogen
Ceruloplasmin	Prealbumin
Clusterin	Retinol Binding Protein
CRP	SAA
Cystatin C	SAP
Ferritin	Thiostatin
Fetuin A	Thymidine Kinase (TK1)
Fetuin	Transferrin
Fibrinogen	Vitamin D Binding Protein
GST alpha	
Hamster (CHO) Annexin A5	
Hamster (CHO) Beta 2-Microglobulin B2M	
Hamster (CHO) C-C Motif Chemokine	
Hamster (CHO) Cathepsin B	
Hamster (CHO) Cathepsin D	
Hamster (CHO) Cathepsin L	
Hamster (CHO) Cathepsin Z	
Hamster (CHO) Clusterin	
Hamster (CHO) GSTp	
Hamster (CHO) GTP-Binding Nuclear Protein (RAN)	
Hamster (CHO) Legumain	
Hamster (CHO) LPLA2	
Hamster (CHO) Lysosomal Acid Lipase	
Hamster (CHO) Lysosomal Protective Protein (CTSA)	
Hamster (CHO) Matrix metalloproteinase-19 (MMP-19)	
Hamster (CHO) Nidogen-1	
Hamster (CHO) Nucleobindin-2	
Hamster (CHO) Peroxiredoxin-1	
Hamster (CHO) PLBL2	
Hamster (CHO) Transforming Growth Factor Beta-1	
Haptoglobin	
Hemoglobin	
Hemopexin	



The ICL Difference

- Robust catalog of antibodies, ELISA kits, and protein reference standards
- Consistent batch-to-batch performance
- Reliable long-term supply
- Expert technical support
- Antigen sequencing reproducibility safeguards
- Rigorous in-house testing and gold standard quality control
- Hundreds of product citations from high-impact scientific journals



Immunology Consultants Laboratory, Inc.
Manufactured in Portland, Oregon, USA

T: 503.747.2454 F: 503.747.2544 WWW.ICLLAB.COM SALES@ICLLAB.COM